


# APESIN spray F



## Alcohol-based rapid disinfection, liquid - perfume-free

- Fast-acting
- Free of dyes and perfume

### Product Profile

- APESIN spray F is a fast-acting disinfectant against bacteria, yeasts, and viruses (limited spectrum).
- This product fulfils the latest hygiene requirements and complies with EU Biocidal Products Regulation (BPR).
- It is suitable for all alcohol resistant surfaces and furniture in healthcare facilities, in industrial and institutional areas as well as in food areas.
- APESIN spray F saves one working step because food contact surfaces do not need to be rinsed with water after the use of the product.
- APESIN spray F is free of chlorine, aldehyde, and phenol.

### Area of application

- The product can be used on mattresses, beds, transport trolleys, instrument deposits, carpets and bathroom tiles, dustbins, toilet seats, frameworks, door handles, furnitures and work surfaces, handrails, shoes, as well as in places that are difficult to access.
- Also suitable for the cleaning of cutter boxes and containers. It is best appropriated for knives, boots and cleaning machines.
- APESIN spray F can be applied on aluminium surfaces.
- Do not apply on alcohol sensitive materials (e.g.: PLEXIGLAS®).
- For painted surfaces, test material resistance on an inconspicuous area first.

### Application note

- Apply undiluted. Spray directly onto objects from a distance of 20-30 cm, so that surfaces are completely wetted, or apply on a clean and absorbent cloth and wipe the surface completely wet.
- Max. 50 ml/m<sup>2</sup>. Do not spray larger surfaces (max. 2m<sup>2</sup>).
- Allow the solution to work, for detailed information and exposure time, see next page: Standards & spectrum of activity.

### Ingredients:

Ingredients per 100g: 40g ISOPROPYL ALCOHOL, 25g PROPYL ALCOHOL

### Sustainable development and production:



### Product safety, storage and environmental instructions

**Safety:** For professional use only. For further information, see Safety Data Sheet. Use biocides safely. Always read the label and product information before use.

**Storage:** Store at room temperature in the original container.

**Environmental:** Only give completely emptied container to special waste disposer.

### Application & Dosing



Follow below instructions.



Conc.

**For disinfection:**  
Suitable for areas with specific hygiene requirements.



Conc.

**Spray application:**  
Spray onto a clean cloth and wipe the surface wet or spray only.



**Surface cleaning:**  
Wipe the surface wet or damp with a clean cloth.



Conc.

**Toilet cleaning:**  
Clean toilet seats.



Conc.

Ideal for application on surfaces in the kitchen area.




### Sales units:

Order No. 714867      10 x 750 ml  
Order No. 714929      10 x 750 ml

### Your trusted partner on site:



# APESIN spray F



## Alcohol-based rapid disinfection, liquid - perfume-free

| Surface Disinfection:<br>Food area / industrial & institutional area | contact time |        |        |        |
|--|--------------|--------|--------|--------|
|  | 30 sec.      | 1 min. | 2 min. | 5 min. |

|                              |                   |  |   |  |
|------------------------------|-------------------|--|---|--|
| <b>bactericidal activity</b> |                   |  |   |  |
| EN 1276                      | high organic load |  | X |  |
| EN 13697                     | high organic load |  | X |  |
| <b>yeasticidal activity</b>  |                   |  |   |  |
| EN 1650                      | high organic load |  | X |  |
| EN 13697                     | high organic load |  | X |  |

| Surface Disinfection:<br>Medical area | contact time |        |        |        |
|---------------------------------------|--------------|--------|--------|--------|
|                                       | 30 sec.      | 1 min. | 2 min. | 5 min. |


|   |                   |   |   |   |
|---|-------------------|---|---|---|
| <b>bactericidal activity</b>                          |                   |   |   |   |
| EN 13727  | high organic load | X |   |   |
| EN 13697  | high organic load |   | X |   |
| EN 16615  | high organic load |   | X |   |
| <b>tuberculocidal activity</b>                        |                   |   |   |   |
| EN 14348  | high organic load | X |   |   |
| <b>yeasticidal activity</b>                           |                   |   |   |   |
| EN 13624  | high organic load | X |   |   |
| EN 13697  | high organic load |   | X |   |
| EN 16615  | high organic load |   | X |   |
| <b>virucidal activity</b>                             |                   |   |   |   |
| EN 14476 Virucidal activity against enveloped viruses | high organic load | X |   |   |
| EN 14476 Limited spectrum virucidal activity          | high organic load |   |   | X |
| <b>EN 14476 individual test organisms</b>             |                   |   |   |   |
| Vaccinia Virus  | high organic load | X |   |   |
| Murine Norovirus (MNV)                                | high organic load |   | X |   |
| Adenovirus  | high organic load |   |   | X |

| VAH-Certificate | contact time |        |        |        |
|-----------------|--------------|--------|--------|--------|
|                 | 30 sec.      | 1 min. | 2 min. | 5 min. |

|   |   |  |  |   |
|---|---|--|--|---|
| Surface disinfection in medical areas, bacterial and yeasticidal activity | with mechanical action dirty conditions |  |  | X |
|---|---|--|--|---|



# APESIN spray F



## Alcohol-based rapid disinfection, liquid - perfume-free

| <b>Surface Disinfection:</b>  |   | <b>contact time</b> |        |        |        |
|---|---|---------------------|--------|--------|--------|
| Food area / industrial & institutional area                               |   | 30 sec.             | 1 min. | 2 min. | 5 min. |
| bactericidal + yeasticidal activity                                       |   |                     |        |        |        |
| EN 1276 + EN 1650 + EN 13697  | high organic load                       |                     | X      |        |        |
| <b>Surface Disinfection:</b>  |   | <b>contact time</b> |        |        |        |
| Medical area  |   | 30 sec.             | 1 min. | 2 min. | 5 min. |
| bactericidal + yeasticidal activity (without mechanical action)           |   |                     |        |        |        |
| EN 13727 + EN 13624+ EN 13697   | high organic load                       |                     | X      |        |        |
| bactericidal + yeasticidal activity (without mechanical action)           |   |                     |        |        |        |
| EN 13727 + EN 13624 + EN 16615  | high organic load                       |                     | X      |        |        |
| tuberculocidal activity   |   |                     |        |        |        |
| EN 14348  | high organic load                       | X                   |        |        |        |
| virucidal activity  |   |                     |        |        |        |
| EN 14476 Virucidal activity against enveloped viruses                     | high organic load                       | X                   |        |        |        |
| EN 14476 Limited spectrum virucidal activity                              | high organic load                       |                     |        |        | X      |
| EN 14476 individual test organisms  |   |                     |        |        |        |
| Vaccinia Virus  | high organic load                       | X                   |        |        |        |
| Murine Norovirus (MNV)  | high organic load                       |                     |        | X      |        |
| Adenovirus  | high organic load                       |                     |        |        | X      |
| <b>VAH-Certificate</b>  |   | <b>contact time</b> |        |        |        |
|   |   | 30 sec.             | 1 min. | 2 min. | 5 min. |
| Surface disinfection in medical areas, bacterial and yeasticidal activity | with mechanical action dirty conditions |                     |        |        | X      |

